

Try out the online ROMA calculator available on the Elecsys® HE4 page at cobas.com

Download the Roche ROMA application for the iPhone® and the iPad® from the App Store.



Roche ROMA



#### References

- 1 Huhtinen, K. et al. (2009). *Br J Cancer*. 100 (8): 1315-1319.
- 2 Moore, R.G. et al. (2012a). *Am J Obstet Gynecol*. 206: 351.e1-8.
- 3 Lu, R. et al. (2012). *Biochem. Biophys. Res. Comm.* 419: 274-280.
- 4 Zhu, Y.-F. et al. (2013). *Asian Pacific Journal of Tropical Medicine*. 6(4): 265-272.
- 5 Moore, R.G. et al. (2012b). *Scientific Reports*. 4, Article number: 3574.
- 6 Paulsen, T. et al. (2006). *Int J Gynecol Cancer*. 16 (Suppl. 1): 11-17.
- 7 Bristow, R.E. et al. (2007). *Cancer*. 109:1513-22.
- 8 Roche HE4 package insert.
- 9 Ortiz-Munoz, B. et al. (2014). *Tumor Biol*. 35(7):7249-58.
- 10 Moore, R.G. et al. (2011). *Obstet Gynecol*. 118 (2, Part 1): 280-288.
- 11 Karlsen, M.A. et al. (2012). *Gynecol Oncol*. 127: 379-383.
- 12 Patient cases kindly provided by Dr. Rafael Molina, President of ISOBM. Oncology Unit, Laboratory of Clinical Biochemistry, Hospital Clinic, Barcelona, Spain.

**The App and the online calculator are for educational use only.**

iPhone and iPad are trademarks of Apple Inc.

COBAS, ELECSYS and LIFE NEEDS ANSWERS are trademarks of Roche.

©2014 Roche

Roche Diagnostics International Ltd  
CH-6343 Rotkreuz  
Switzerland  
www.cobas.com

**Assessing ovarian cancer risk with the ROMA algorithm**  
*Combining CA 125 and HE4 tumor markers in pelvic mass evaluation*




**cobas**<sup>®</sup>  
Life needs answers

# Introduction of HE4 tumor marker


*HE4 level is highly expressed in epithelial ovarian cancer*

In a study by Huhtinen, K. et al. (2009), serum CA 125 and HE4 concentrations of healthy patients were compared with patients with endometriosis (Endo) and epithelial ovarian cancer (EOC)<sup>1</sup>

- HE4 level is highly elevated in epithelial ovarian cancer
- HE4 shows high differentiation between ovarian cancer vs. endometriosis and healthy controls




*HE4 is involved in ovarian cancer cell proliferation and invasion<sup>3-5</sup>*




*Patients with high HE4 levels have poorer outcomes, suggesting HE4 plays an important role in ovarian cancer*

In a study by Moore, R.G. et al. (2012a), serum were obtained from 1,042 women with benign diseases prior to surgery for pelvic mass and HE4 levels were determined.<sup>2</sup>

- HE4 levels were in general not elevated in benign diseases, including ovarian cysts and endometriosis



**Patients with Stage I-IV EOC (N=89)**




In a study by Moore, R.G. et al. (2012b), Kaplan-Meier survival curves of 89 ovarian cancer patients (all stages) were shown.<sup>5</sup>

The graph shows decreased survival for patients with HE4  $\geq 500$  pmol/L compared to patients with HE4  $< 500$  pmol/L. Women with high HE4 levels had a 5-year median overall survival of 27% compared to 59% for those with low HE4 levels. Hazard ratio = 2.2 ( $p=0.005$ ).

# Accurate pelvic mass evaluation leads to accurate referral decision


## Referral decision impacts patient survival and cost effectiveness of treatments

Survival rate in ovarian cancer patients is dependent on the type and the quality of surgery to remove the tumor as completely as possible. Paulsen et al. (2006) showed that, among 198 women with diagnosis of advanced epithelial invasive ovarian, tubal, and peritoneal cancer, the highest survival rate can be achieved when treatment is done with gynecological oncologists and in expert centers (e.g. teaching hospitals). This improved survival may be because the surgery for ovarian cancer is complicated and gynecological oncologists or specialists at expert centers are trained to perform the comprehensive debulking and surgical staging needed for optimal treatment of ovarian cancer.<sup>6</sup>



Women who had surgery performed by gynecological oncologists had 20% increase in short-term survival compared to women who had surgery performed by non-gynecological oncologists ( $p < 0.0001$ ).<sup>6</sup>

Bristow et al. (2007) showed that referral of patients with ovarian cancer to expert centers was more cost-effective than referral of patients with ovarian cancer to less experienced centers: \$9,893 per quality-adjusted life year (QALY) compared with \$17,149 per QALY, respectively.<sup>7</sup>



Women who were operated at teaching hospitals had longer survival compared to those who were operated at non-teaching hospitals (Hazard ratio=1.83;  $p=0.02$ ).<sup>6</sup>


# ROMA, Risk of Ovarian Malignancy Algorithm, in pelvic mass evaluation

## CA 125 and HE4 can be combined in a mathematical algorithm to better assess the risk of epithelial ovarian cancer in women with pelvic mass

### Risk assessment with the ROMA (Risk of Ovarian Malignancy Algorithm) score:

The score is used to calculate the risk of having epithelial ovarian cancer in women arriving at the physician with pelvic mass. It includes the biomarker levels of HE4 and CA125 and also considers the menopausal status of the woman.

The dual marker combination can be used to classify women into high and low risk groups, allowing for the effective triage of women to appropriate surgical centers for their care.<sup>2,8</sup>




Calculation of the ROMA scores for pre- and post-menopausal women and individual cut-points for the Roche assays to separate between low and high risk patients<sup>9</sup>

# Clinical benefits of ROMA


*Compared to testing CA 125 alone, ROMA has higher specificity for ovarian cancer detection*

Ovarian cancer vs. benign gynecologic diseases



This means ...


**Among 10 women with benign gynecologic diseases, ROMA dismisses 1 more patient which CA125 alone would include as having ovarian cancer**



*ROMA accurately identifies 94% (45/48) of the patients with pelvic mass that have ovarian cancer*

In a study by Moore, R.G. et al. (2011), a total of 472 patients were evaluated: 383 women with benign diseases and 48 women with epithelial ovarian cancers (9 pre-menopausal and 39 post-menopausal)<sup>10</sup>


- Use of HE4 and CA 125 with ROMA has a high sensitivity for the prediction of epithelial ovarian cancer in women with pelvic mass



EOC = Epithelial ovarian cancer


*Compared to testing CA 125 alone, ROMA has higher sensitivity for ovarian cancer detection*

Ovarian cancer vs. benign gynecologic diseases



This means ...

**Among 10 women with ovarian cancer, ROMA identifies 1 more patient than CA 125 alone**




*Compared to testing CA 125 alone, ROMA increases sensitivity by 10% in the detection of stages I/II ovarian cancer*

In a study by Karlsen, M.A. et al. (2012), serum from 1,218 patients in the prospective pelvic mass study was collected prior to diagnosis.<sup>11</sup>

- Diagnosis: 809 benign tumors, 79 borderline ovarian tumors, 252 ovarian cancer (64 early and 188 late stage), 9 non-epithelial ovarian tumors and 69 non-ovarian tumors
- The study concluded that ROMA increases differentiation of ovarian cancer from other pelvic masses than CA 125 alone, even in early stage ovarian cancer
- ROMA might be valuable as a first line biomarker for selecting high risk patients for referral to a tertiary center and further diagnosis

Stages I/II ovarian cancer vs. benign gynecologic diseases



In a study by Ortiz-Munoz, B. et al. (2014), pre-operative serum levels of HE4 and CA 125 were measured in 119 women with benign gynecologic diseases, 67 patients with ovarian cancer and 32 healthy women.<sup>9</sup>

- ROMA has higher specificity and sensitivity than CA 125 alone to discriminate ovarian cancer from benign gynecologic diseases.

## ROMA in action: Patient cases<sup>12</sup>

### Patient case #1

A 44-year-old (pre-menopausal) female with abdominal pain and a history of endometriosis was detected by ultrasound to have a pelvic mass. Her test results were the following:

Parameter	Result	Normal interval	Report: High CA 125 level but normal HE4 level Diagnosis: Endometrial cyst and exacerbation of endometriosis
Creatinine	0.7	<0.9 mg/dL	<p style="text-align: center;"><b>ROMA score: Low risk</b></p> <p style="text-align: center;">Low risk ROMA score supports the diagnosis of non-malignant tumor</p>
AST	22	<35 U/L	
ALT	19	<35 U/L	
CA 19-9	22	<39 U/mL	
CA 125	150	<35 U/mL	
HE4	34	<140 pmol/L	
ROMA	3.6	<11.4% for pre-menopausal	

### Patient case #3

A 58-year-old (post-menopausal) female, non-smoker with no toxic habits was referred following the detection of a pelvic mass during an ultrasound scan. Her test results were the following:

Parameter	Result	Normal interval	Report: High CA 125 and HE4 levels Diagnosis: Ovarian cancer
Creatinine	0.86	<0.9 mg/dL	<p style="text-align: center;"><b>ROMA score: High risk</b></p> <p style="text-align: center;">High risk ROMA score supports the diagnosis of ovarian cancer</p>
AST	32	<35 U/L	
ALT	23	<35 U/L	
CA 19-9	6	<39 U/mL	
CA 125	324	<35 U/mL	
HE4	467	<140 pmol/L	
ROMA	92.6	<29.9% for post-menopausal	

### Patient case #2

A 52-year-old (pre-menopausal) female was detected by ultrasound to have a pelvic mass. Her test results were the following:

Parameter	Result	Normal interval	Report: High HE4 level but normal CA 125 level Diagnosis: Ovarian cancer
Creatinine	0.7	<0.9 mg/dL	<p style="text-align: center;"><b>ROMA score: High risk</b></p> <p style="text-align: center;">High risk ROMA score supports the diagnosis of ovarian cancer</p>
AST	22	<35 U/L	
ALT	45	<35 U/L	
CA 19-9	6	<39 U/mL	
CA 125	19	<35 U/mL	
HE4	601	<140 pmol/L	
ROMA	96.8	<11.4% for pre-menopausal	

### Patient case #4

A 59-year-old (post-menopausal) female, non-smoker was found to have a pelvic mass. Her test results were the following:

Parameter	Result	Normal interval	Report: Normal CA 125 and HE4 levels Diagnosis: Ovarian cyst
Creatinine	0.9	<0.9 mg/dL	<p style="text-align: center;"><b>ROMA score: Low risk</b></p> <p style="text-align: center;">Low risk ROMA score supports the diagnosis of non-malignant tumor</p>
AST	32	<35 U/L	
ALT	25	<35 U/L	
CA 19-9	16	<39 U/mL	
CA 125	29	<35 U/mL	
HE4	83	<140 pmol/L	
ROMA	26.3	<29.9% for post-menopausal	

